

Thailand Bioresource Research Center (TBRC)

National Center for Genetic Engineering and Biotechnology

Innovation Cluster 2 (Tower B, 8th Floor)

143 Thailand Science Park, Phahonyothin Road

Khlong Nueng, Khlong Luang, Pathum Thani 12120, Thailand

Instructions for sending microorganisms for preservation service at the Thailand Bioresource Research Center (TBRC)

1. Preparation of documents

The service form must be submitted to TBRC together with the microorganisms.

2. Preparation of microorganisms

- 2.1. Grow microorganisms on suitable/sporulation dry-surface agar media to reduce risk of contamination during shipment.
 Send 1 tube or plate per strain of microorganism.
 - 2.1.1. Fast growing microorganisms (most bacteria and yeasts) should be sent within 24 hours after inoculation.
 - 2.1.2. Sporulating filamentous fungi and actinomycetes should be sent after spores are mature. Non-sporulating fungi (can be preserved using freezing method only) should be sent within 1-2 weeks after inoculation.
 - 2.1.3. The customer can send microorganisms in liquid-dried or freeze-dried form. Service fee will be applied for cultivation of microorganisms from liquid-dried or freeze-dried ampoules. Send 1 tube per strain of microorganism.
 - 2.1.4. The customer may be requested to send culture media with microorganisms.

3. Shipment of microorganisms to TBRC

3.1. Drop-off at TBRC

Please contact TBRC in advance (Tel: +66 (0)2 117 8001, email: tbrcservice@biotec.or.th). Microorganisms should be brought to the TBRC office located at Innovation Cluster 2 Tower B, 8th Floor, Thailand Science Park, Phahonyothin Road, Khlong Nueng, Khlong Luang, Pathum Thani.

3.2. By post

- 3.2.1. Microorganism tubes/plates must be packed in plastic bags or plastic tubes and properly sealed before putting in a cardboard box for shipment (do not pack in an envelope or a paper bag). Microorganisms classified in Risk Group 2 must be packed and shipped according to national biosafety and biosecurity regulations or IATA Dangerous Goods Regulations.
- 3.2.2. Send microorganisms, service form and proof of customer identity by express mail service (EMS) to the address indicated below.

То

Preservation Service
National Center for Genetic Engineering and Biotechnology
BIOTEC Building, 3rd Floor, Room 313
113 Thailand Science Park, Phahonyothin Road
Khlong Nueng, Khlong Luang, Pathum Thani 12120, Thailand